



# HOT PEEL TRANSFERS

**\* SAMPLE TRANSFER INCLUDED-CUT UP & TEST FIRST \***

**1. THREAD GARMENT ON LOWER PLATEN**

**2. PREHEAT GARMENT**

**3. PRESS TRANSFER**

\*INSTRUCTIONS BASED ON OUR HEAT PRESSES MADE AT FIRST EDITION.  
YOU MAY NEED TO ADJUST YOUR TIME AND/OR TEMP BASED ON YOUR HEAT PRESS.

 **HEAVY (50-60 PSI)**

 **380°**

 **15-17 SECONDS**

NOTE - MAY NEED TO INCREASE TEMP/ TIME TO CURE TRANSFER

**4. DELAY PEEL FOR 3-5 SECONDS**